

MEP Instruments Media Release

“New automation possibilities in Ion Chromatography”

For immediate publication

Number of words: 205

Media & Technical Contact:

Bart Cleeren, MEP Instruments Pty Ltd, bc@mep.net.au

Photo: Metrohm 815 Robotic Soliprep for IC

Keywords: Ion Chromatography, Automation, Sample Preparation, Liquid Handling

Metrohm

New automation possibilities in ion chromatography: Metrohm 815 Robotic Soliprep

Pressure to increase output while reducing costs requires laboratory automation wherever possible. Metrohm’s new 815 Robotic Soliprep meets such demands as it allows to completely automate preparation and analysis of liquid and solid samples in ion chromatography. The 815 Robotic Soliprep provides maximum flexibility as it can be customized on demand to meet individual customers’ specific requirements.



The 815 Robotic Soliprep is a state-of-the-art sample processor, which reliably transfers both liquid and solid samples to an ion chromatograph. A broad range of working steps such as pulverizing, homogenizing, extracting, filtering or diluting can be safely entrusted to the 815 Robotic Soliprep and included in the analysis method. Due to the incorporation of inline sample preparation into the analysis system, manual intervention is no longer required during the determination. The sample rack of the 815 Robotic Soliprep is equipped with various inserts, sample beakers, syringe filters and needles and can be customized on demand to meet individual customers’ specifications. Combined with Metrohm’s advanced ICs the 815 Robotic Soliprep opens up numerous

new applications and application fields. Applications include, to mention just a few, the analysis of **ingredients and contaminants in tablets, animal feed, sediments or food.**

Read more about Metrohm 815 Soliprep at www.metrohm.com

**MEP Instruments Pty Ltd, PO Box 1880, North Ryde 2113,
Phone 02 98786900. Fax: 02 98786999**

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Analytik Jena, Applikon, Autolab, Camag, Retsch and Retsch Technology in Australia and New Zealand.

Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113
Tel 02 9878 6900 Fax 02 9878 6999 info@mep.net.au; or in New Zealand:

MEP Instruments Ltd (New Zealand) 2 B, Unit D, Building 2, 5 Ceres Court, North Shore City 0630
PO Box 305161, North Shore 0757
Phone: (+64) 09 912 1330 Fax: (+64) 09 912 6892 info@mep.net.nz Visit MEP online:
<http://www.mep.net.au>

