

## Press Release

for immediate publication

### Ecotrode Gel – maintenance-free gel electrode with aging indicator

**Metrohm introduces the Ecotrode Gel – the new maintenance-free electrode for pH measurement and titration. The new aging indicator tells users at a glance, if the sensor is still good for use. The main properties of the Ecotrode Gel:**

- maintenance-free electrolytes: refilling or checking the electrolyte not necessary
- aging indicator: decolorization of the white electrode gel tells users exactly when replacement of the sensor is due
- ideal for routine measurements



Image: 2843\_s.jpg

Weblink: [news.metrohm.com](https://news.metrohm.com)

Keywords: pH measurement, titration, electrode maintenance

Branches: all kinds of industry

**MEP Instruments Pty Ltd, PO Box 1880, North Ryde 2113,  
Phone 02 98786900. Fax: 02 98786999**

MEP Instruments is sole distributor of leading names in sample analysis technology in Australia and New Zealand. The joint venture company formed by Metrohm Switzerland and Anton Paar Austria provides comprehensive service, support and training for the whole product line of Metrohm, Anton Paar, Analytik Jena, Applikon, Autolab, Camag, Retsch and Retsch Technology in Australia and New Zealand.

Contact MEP Instruments Pty Ltd (Australia), Unit 5, 112-118 Talavera Road, North Ryde, NSW 2113 Tel 02 9878 6900 Fax 02 9878 6999 [info@mep.net.au](mailto:info@mep.net.au); or in New Zealand:

MEP Instruments Ltd (New Zealand) 2 B, Unit D, Building 2, 5 Ceres Court, North Shore City 0630

PO Box 305161, North Shore 0757

Phone: (+64) 09 912 1330 Fax: (+64) 09 912 6892 [info@mep.net.nz](mailto:info@mep.net.nz) Visit MEP online: <http://www.mep.net.au>